

a beautiful me

AO+ SKINCARE VITA PEN

AO+ Vita Pen is a cordless and rechargeable device that uses the 630nm Low Power Laser Light (LPLL) and low current electric pulse to relieve pain and accelerate the healing of skin condition. It also helps to repair and regenerate cells and damaged tissues.

Beauty

Use Vita-Pen with AO+ Geranium Ozonated Olive Oil to improve blood circulation, reduce wrinkles, fine lines, as well as to boost collagen production for healthier-looking skin. You can now enjoy the daily facial treatment at the convenience of your home.

(Please see diagrams for facial treatment)

Rashes & Eczema

The use of Vita-Pen with AO+ Ozonated Olive Oil helps to accelerate natural healing for Eczema skin conditions.

Headache/ Migraine

Apply AO+ OZO-10 on the temple and the back of the neck, then use Vita Pen to massage over it.


AO+

Disclaimer: Beware of unauthentic crystal roller made of dyed marble may release toxins and cause harm to the skin.

a beautiful me

Directions of use

- Press and hold the power button until you hear a beep
- Adjust the electric-pulse strength from 1 to 9, to your preference
- Apply AO+ Ozonated Olive Oil on the skin and use it to massage over it
 - Clean the probe after use
- Use the cable provided for recharging the battery
 - **Dome Head - for beauty effects*
 - **Spheroidal Head - for pain relief*



Reduces wrinkles and tightens skin

Stimulates the production of collagen and elastin

Reduces inflammation and signs of aging

Promotes healthy blood circulation

Speeds up healing processes for rashes

Relieves itch & muscle pain

Provides superficial therapy & relieves joints pain effectively

Quick relief for migraine

Promotes hair growth



aoplus.sg

AO+